
Deutsch	4
English	58

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

for purchasing this Ondal product. Please read these Operating Instructions very carefully, abide by the safety notices and observe all operating and cleaning requirements.

For which appliances do these Operating Instructions apply?

- ACROBAT 77/2000 pendant system equipped with:
 - **spring arm (standard)** 
 - **spring arm (NRH)** 
- ACROBAT LCH equipped with **spring arm (LCH)** :
 - Ceiling mounted version:
 - Single type – double-arm,
 - Dual type – double-arm,
 - Wall mounted version:
 - Type – two arms,

Please do not hesitate to contact our Customer Service team

if you have any questions about the appliance and its installation, and also in service or warranty cases.

Manufacturer

Ondal Medical Systems GmbH • Wellastraße 6 • D-36088 Hünfeld • Germany
Phone: +49 / (0)6652 / 81-600
Fax: +49 / (0)6652 / 81-392
© Ondal Medical Systems GmbH, 2012
1505093, Edition 2012-05, Version 5


Visit us on the Internet

www.ondal.de

E-mail address

info@ondal.de

Supplier's address

Space for supplier's stamp or label

Notes on copyright and property rights

- Copyright
- All rights reserved. These Operating Instructions are protected by copyright law.
 - Any other use than that regulated by law must be approved in writing by Ondal Medical Systems GmbH.
 - Ondal does not take any liability for or in relation to the misuse of this information in the prohibited manner by any person or company.

Modifications and translations

- Modifications to the appliance
- Ondal products are subject to continuous further development. Ondal reserves the right to modify the form, equipment and technology of the scope of supply without prior notice.
- Modifications to the Operating Instructions
- The contents of these Operating Instructions are subject to change without prior notice.
- Translations
- In case of translations into foreign languages, the German version of these Operating Instructions shall take precedence.

Trademarks

- All trademarks mentioned in these Operating Instructions are the sole and exclusive property of the corresponding manufacturer.

1	Important Information for Safe Use	62
1.1	Information for identification of the device	62
1.2	How to identify the Operating Instructions	62
1.3	Identification of target groups	62
1.3.1	Operator	62
1.3.2	User	62
1.3.3	Qualified personnel	62
1.4	Notes for the operator	63
1.4.1	Initial commissioning	63
1.4.2	Availability of these Operating Instructions	63
1.4.3	Warranty	63
1.5	Notes for the user	64
1.5.1	Instruction on the device	64
1.5.2	User's duty to inform and inspect	64
1.5.3	Marking	64
1.5.4	Standards and guidelines	64
1.6	Intended purpose	65
1.6.1	Incorrect use	65
1.6.2	Contraindications	65
1.7	Ambient conditions for operation and storage	65
1.7.1	Ambient conditions for storage and transport	65
1.7.2	Ambient conditions for operation	65
1.8	Approved Ondal adaptations	66
1.9	Combination with products of other manufacturers	66
1.10	Patient environment	66
1.11	Disposal	67
2	Safety Instructions	68
2.1	Structure of the safety instructions	68
2.1.1	Warnings of risk of injury	68
2.1.2	Warnings of damage to property	68
2.1.3	Indication of additional information	68
2.2	Supplementary symbols used in the safety instructions	68
2.3	Graphic symbols on the appliance and/or on the packaging	69
2.4	Information on the rating plate	70
2.5	Overview of the most important safety instructions	70
2.5.1	Operation	70
2.5.2	Mounting / dismantling	71
2.5.3	Cleaning and disinfection	71
2.5.4	Maintenance work	72
3	Device and Functional Description	73
3.1	Parts and control elements	73
3.2	Validity of these Operating Instructions	74
3.3	Appliance description	74
3.4	Functional description	75
4	Initial and Repeated Commissioning	76
5	Operation	77
5.1	Positioning the pendant system	77
5.1.1	Swiveling the pendant system	77
5.1.2	Adjusting the height of the pendant system	77
6	Cleaning and Disinfection	78
6.1	General safety instructions	78
6.2	Cleaning	78
6.3	Disinfection	78

7	Maintenance	79
7.1	General safety instructions	79
7.2	Built-in components from third-party manufacturers	79
7.3	Maintenance of the securing segment	79
7.4	Repeated inspection	79
8	Dismantling / Mounting in Case of Service	80
8.1	General safety instructions	80
8.2	Calculating the maximum loading capacity (payload)	81
8.2.1	What is the maximum loading capacity (payload)?	81
8.2.2	Determining the maximum loading capacity (payload)	81
8.2.3	Calculating the maximum loading capacity (payload) of the ACROBAT 2000 spring arm	82
8.2.4	Calculating the maximum loading capacity (payload) of the NRH spring arm	83
8.2.5	Calculating the maximum loading capacity (payload) of the LCH spring arm	84
8.3	Dismantling the end device from the ACROBAT2000 spring arm	85
8.4	Mounting the end device to the ACROBAT2000 spring arm	86
8.5	Dismantling the end device from the NRH spring arm	88
8.6	Mounting the end device to the NRH spring arm	89
8.7	Dismantling the end device from the LCH spring arm	91
8.8	Mounting the end device to the LCH spring arm	92
9	Adjustments	95
9.1	General safety instructions	95
9.2	Adjusting the spring tension	95
9.2.1	Adjusting the spring tension on the ACROBAT2000 and NRH spring arm	95
9.2.2	Adjusting the spring tension on the ACROBAT LCH spring arm	96
9.3	Adjusting the brake force	97
9.3.1	Adjusting the brake force on the pendant system	97
9.3.2	Adjusting the brake force on the ACROBAT2000 NRH spring arm	98
9.3.3	Adjusting the brake force on the ACROBAT LCH spring arm	99
9.4	Adjusting the vertical lift	100
9.4.1	Adjusting the vertical lift on the ACROBAT 2000 spring arm	100
9.4.2	Adjusting the vertical lift on the ACROBAT 2000 NRH spring arm	100
9.4.3	Adjusting the vertical lift on the ACROBAT LCH spring arm	101
10	Technical Description	102
10.1	ACROBAT DB 2001 ceiling-mounted version, single-arm variant	102
10.2	ACROBAT DB 2002 ceiling-mounted version, dual-arm variant	103
10.3	ACROBAT wall-mounted version, single-arm variant	104
11	Technical Data	105
12	Optional Accessories	107
13	Inspection Plan	108

Device identification

1.1 Information for identification of the device

- These Operating Instructions are intended solely for devices with the manufacturer's rating plate bearing the following information:
 - Device designation: ACROBAT 2000 pendant system
 - Type designation: ACROBAT 2000, ACROBAT LCH

Make sure you are using the latest version

Identification of these Operating Instructions

1.2 How to identify the Operating Instructions

- To ensure that you always have the latest version of these Operating Instructions, all pages bear a 7-digit identity number including the date of issue and the version number:
 - Edition: 5/8/2012
- This identification is binding for the validity of the Operating Instructions and must not be removed, regardless of the type of publication (printed form, electronic form or excerpts).

1.3 Identification of target groups

The groups of persons described below are mentioned in these Operating Instructions:

1.3.1 Operator

The following natural persons or legal entities shall be considered as operators:

- All persons who use the appliance in a medical practice, hospital, etc. or hand over the appliance to third parties for use/application, and who have actual physical authority over the appliance during operation.
- The operator shall be liable for handing over a safe appliance and for instructing the user in its proper operation and normal use.

1.3.2 User

The following persons shall be considered as users:

- Persons who, due to their professional qualification and instruction by the persons designated by the operator, are authorised to operate the appliance and to work with it.
- Users shall be fully responsible for the safe operation of the appliance in accordance with its intended purpose.

1.3.3 Qualified personnel

The following persons shall be considered as qualified personnel:

- Persons who underwent special professional training in the field of medicine or medical engineering,
- persons who can assess their work and recognise the potential hazards involved on the basis of their professional experience and instruction in safety-relevant regulations.
- In States where the performance of tasks in the medical or medical engineering sector are subject to certification, qualified personnel must have obtained the corresponding certificate.

1.4 Notes for the operator

- Although the appliance has been designed according to the state of the art and is safe to operate, the appliance must be considered a potential source of danger, in particular when operated by insufficiently trained personnel or used improperly and not as prescribed.
- The appliance may only be operated, cleaned and disinfected by trained qualified personnel.
- All the mounting, dismantling and adjustment work described in these Operating Instructions may only be carried out by qualified personnel which have been authorised and instructed by the operator.
- For safety reasons, any operator actions or interventions exceeding this scope may only be carried out by Ondal or companies authorised by Ondal. As a prerequisite for the authorisation of a company, its service technicians must have successfully participated in technical training organised by Ondal. This authorisation is granted for a limited period.

1.4.1 Initial commissioning

Validity

- These Operating Instructions only apply after initial commissioning has been carried out properly.
- Prior to initial use, the appliance must be thoroughly cleaned and disinfected.
- The instruction for the proper installation of the appliance is included in the Installation Instructions applicable for the appliance.

1.4.2 Availability of these Operating Instructions

Duty to inform

- Since these Operating Instructions are an integral part of the appliance, they must always be kept near the appliance in order to be able to look up safety instructions and important information on use at any time.
- Do not pass on the appliance to any third party without valid Operating Instructions. Based on the ID and version numbers, make sure that you hand over an up-to-date and valid version of the Operating Instructions together with the appliance.

1.4.3 Warranty

The warranty of Ondal for the safety and operational reliability of the appliance is subject to the following conditions:

- The appliance is used exclusively as prescribed and operated as specified in these Operating Instructions.
- Only genuine spare parts or accessories and those defined and approved by Ondal Medical Systems GmbH are used. The use of other parts may involve unknown risks and must be avoided in all cases.
- No structural alterations are made to the appliance. Unauthorised modifications or conversions of the appliance are not permitted for safety reasons.
- Inspections and maintenance must be carried out at the specified time intervals.
- Initial commissioning has been carried out, and the appliance has been released for operation by means of a declaration of acceptance.

1.5 Notes for the user

- All the work steps described in these Operating Instructions may only be carried out by qualified personnel which have been authorised and instructed by the operator.

1.5.1 Instruction on the device

Instruction

- The instruction must be carried out immediately and directly on the appliance by Ondal, by a company authorised by Ondal or by a person designated by the operator.
- On completion of the instruction, a certificate must be created and signed in order to document that the user has understood the special operator control actions required for normal use.

1.5.2 User's duty to inform and inspect

Duty to inform and inspect

- Read these Operating Instructions carefully prior to initial use of the appliance. This ensures that you benefit from all the advantages of the appliance and prevents any risk of injury or damage to property.
- Prior to any use or transfer for use, the functional reliability and proper condition of the appliance must be inspected by the user.

Special problems

- In case of special problems which are not sufficiently described in detail in these Operating Instructions, contact your supplier for your own safety.

1.5.3 Marking

CE mark: Ondal declares that the products comply with the relevant regulations set forth in the applicable European Directives.

Recognised cUL component: This component has been recognised by Underwriters Laboratories Inc. Representative samples of this product have been reviewed by UL and comply with the applicable requirements.

1.5.4 Standards and guidelines

The appliance complies with the safety requirements of the following standards and directives:

- Medical Devices Act (Medizinproduktegesetz - MPG);
- MDD 93/42/EEC, 2007 – Medical Device Directive;
- EN 60601-1: 2006 - Medical Electrical Equipment - Part 1: General Requirements for Basic Safety and Essential Performance.

1.6 Intended purpose

Qualified personnel

- The ACROBAT 2000 pendant system is individually equipped; depending on the version and equipment, it serves for:
 - carrying and positioning medial appliances on the **adaptions** [7] - [9] in the OR, in intensive care and in the medical environment,
- the conveying of media:
 - electricity and data.
- The appliance is suitable for continuous operation.
- The pendant system may only be operated by instructed, qualified medical personnel.
- The pendant system may only be cleaned and disinfected by instructed hygiene specialists.
- Maintenance work on the pendant system must be carried out by the technical specialist personnel of the operator in accordance with the applicable instruction document.

1.6.1 Incorrect use

- The maximum loading capacity of the appliance and its components as specified in *“Chapter 11, Technical Data, on page 105”* must not be exceeded.
- Independently of the maximum loading capacity (payload) the optional components may only be loaded with the maximum weight specified in *“Chapter 12, Optional Accessories, on page 107”*.

1.6.2 Contraindications

- The ACROBAT 2000 pendant system must not be used close to strong magnetic fields.
- No BF or CF application parts in accordance with IEC 60601-1 may be directly connected to the pendant system.

1.7 Ambient conditions for operation and storage

Different conditions are required for the operation and interim storage of the appliance.

1.7.1 Ambient conditions for storage and transport

The following storage conditions apply for storage times of up to 15 weeks:

- Ambient temperature: -25 °C to 70 °C;
- Relative humidity: 10 % to 75 %;
- Atmospheric pressure: 500 hPa to 1060 hPa.

Store only in indoor rooms; after this time, the values specified for the ambient conditions for operation apply.

1.7.2 Ambient conditions for operation

- Ambient temperature: 10 °C to 40 °C;
- Relative humidity: 30 % to 75 %;
- Atmospheric pressure: 700 hPa to 1060 hPa.

Approved accessories

Observe the Operating Instructions for combined products

1.8 Approved Ondal adaptations

- The following Ondal products are approved as accessories to the pendant system:
- Ondal accessories as described in **“Chapter 12, Optional Accessories, on page 107”**.
 - The components are adapted to each other and safe to operate. Any other type of installation, and in particular the use of components from third-party manufacturers, is strictly prohibited because these components can be potential sources of danger.
 - The combination of any other Ondal product with the pendant system must be approved by Ondal Medical Systems GmbH. If applicable, the conformity evaluation must be renewed.

1.9 Combination with products of other manufacturers

- The pendant system can be combined with the adaptations and end devices (e.g. flat screen, OR lamp, etc.) of third-party manufacturers. To prevent dangerous overload which can damage or lead to damage or collapse of the pendant system, the maximum load bearing capacity specified in **“Chapter 11, Technical Data, on page 105”** must be adhered to.
- The party placing the appliance into operation is responsible for the validation of the overall system. A conformity evaluation procedure shall be executed if required and a declaration in accordance with Article 12 of 93/42/EEC (Medical Device Directive, MDD) shall be provided.
- Read the Operating Instructions provided by the third-party manufacturer and in particular the relevant pages with information on the operation of the end device.
- Power packs intended for the supply of end devices must ensure electrical isolation and provide 2 protective measures in accordance with EN 60601-1.


1.10 Patient environment

(See “Figure 1”)

Medical electrical devices which contain tangible parts and are attached to the spring arm and positioned within the patient environment must provide two Means of Patient Protection (MOPP) in accordance with EN 60601-1. If tangible parts are positioned outside the patient environment, two Means Of Operator Protection (MOOP) in accordance with EN 60601-1 must be provided.

The dimensions in the figure show minimum extent of the patient environment in a free surrounding.

Abbildung 1: Patient environment, Figure A.9 from IEC 60601-1:2005

1.11 Disposal

WARNING

Sudden release of spring tension

The spring arm is equipped with a compression spring. When dismantling the spring arm, the compression spring suddenly releases its tension and can lead to severe injury:

- Do not dismantle the spring arm for disposal.

RoHS conformity

- The appliance complies with the requirements of the EU Directive 2002/95/EC ("RoHS Directive" on the restricted use of certain hazardous substances in electrical and electronic appliances).
- To prevent environmental damage and personal injury, we therefore request you to contact us or your authorised service partner if you intend to take the appliance out of operation for the purpose of disposal.
- The appliance must be disposed of at a suitable collection point for the recycling of electrical and electronic devices in accordance with country-specific regulations.

 **DANGER**
 **WARNING**
 **CAUTION**
NOTICE**NOTE**

2.1 Structure of the safety instructions

2.1.1 Warnings of risk of injury

Important notes in these Operating Instructions are marked with graphic symbols and signal words.

Signal words such as DANGER, WARNING or CAUTION describe the degree of risk of injury. The different triangle symbols visually emphasise the degree of hazard.

DANGER refers to an immediate hazard, which, if not avoided, will lead to death or severe injury.

WARNING refers to a potentially dangerous situation which, if not avoided, can lead to death or severe injury.

CAUTION refers to a potential hazard, which, if not avoided, can lead to minor or slight injury.

2.1.2 Warnings of damage to property

NOTICE refers to a potential hazard, which, if not avoided, leads to damage to property.

2.1.3 Indication of additional information

A **NOTE** provides additional information and useful tips for the safe and efficient use of the appliance.

2.2 Supplementary symbols used in the safety instructions

Electric shock hazard: warns of electric shock which can lead to severe injury or even death.

Pendant system crashing: warns of the risk of the pendant system suddenly crashing because the maximum loading capacity has been exceeded.

Sudden release of spring arm: warns that the spring arm may jump up suddenly while dismantling the end device.

2.3 Graphic symbols on the appliance and/or on the packaging

Observe the Operating Instructions: Read these Operating Instructions carefully prior to initial operation of the pendant system. This ensures that you benefit from all the advantages of the pendant system and prevents any risk of injury or damage.

Observe the maximum loading capacity: warns of exceeding the maximum loading capacity (payload) on the pendant system of the **adaption** [7] - [9] or the end device (e.g. flat screen, OR lamp, etc.).

kg: warns of the risk of the **spring arm** [4] suddenly crashing because the maximum loading capacity has been exceeded.


Max. load Nm: warns of the **NRH spring arm** [5] suddenly crashing because the maximum loading capacity has been exceeded.


CE mark: Ondal declares that the products comply with the relevant regulations set forth in the applicable European Directives.


Recognised cUL component: This component has been recognised by Underwriters Laboratories Inc. Representative samples of this product have been reviewed by UL and comply with the applicable requirements.

Atmospheric pressure: indicates the permissible atmospheric pressure values in a range from 500hPa to 1060hPa for transport and storage.

Relative humidity: indicates the permissible humidity values in a range from 10% to 75% for transport and storage.

Ambient temperature: indicates the permissible ambient temperature values in a range from -25°C to 70°C for transport and storage.

Figure 2: Information on the rating plate

TYPE: ACROBAT 2000
 AC 100V 12A 50/60Hz
 AC 120V 12A 50Hz
 AC 230V 13A 50/60Hz

REF 1899483
 SN 1411H000000000

 max. load 21 kg

2011  C  US
 E212032

 Country of Origin Germany
 Ondal Medical Systems GmbH, 36088 Hünfeld

TYPE: ACROBAT LCH
 230V, 13A, 50-60Hz
 120V, 12A, 60 HZ
 24V, 15A

REF 1931833
 SN 1411H000000000

 max. load 176 Nm

2011  C  US
 E212032

 Country of Origin Germany
 Ondal Medical Systems GmbH, 36088 Hünfeld

2.4 Information on the rating plate

The rating plate is attached to the top side of the extension arm (towards the wall bearing or ceiling mount fixture). For the AC 2000 spring arm: under the cover of the supporting tube and on the outside of the cover. For the LCH spring arm: on the top side of the tube profile and under the rear cover.

The illustration shows examples for:

- the ACROBAT 2000 spring arm
- the ACROBAT LCH spring arm

Serial number

- The rating plate indicates the serial number (SN) of the spring arm.

Power supply of the extension arm

- The rating plate provides information on the power supply of the individual spring arms.

Load bearing capacity of the spring arms

- For example, 21 kg or 176 Nm indicates the maximum permissible load bearing capacity of the spring arm/safe working load (see *"Chapter 8.2" on page 81*).

Date of manufacture

- The digits 1 to 4 of the serial number (SN) indicate the date of manufacture of the pendant system.
 - The first two digits indicate the week of manufacture, e.g. 14 = calendar week 14.
 - The following two digits indicate the year of manufacture, e.g. 11 = 2011.
 - The letter in the 5th position indicates the factory, e.g. H = Hünfeld.
 - The digits following the letter indicate the serial number.
- The information and illustrations serve as examples.

2.5 Overview of the most important safety instructions

The safety instructions in the following chapters must be adhered to.

2.5.1 Operation

WARNING

Risk of pendant system crashing because the maximum loading capacity (payload) has been exceeded

If the maximum permissible loading capacity (payload) is exceeded, the pendant system or components of the pendant system may disengage from the fastening device and crash:

- The maximum loading capacity (payload) of the pendant system and components thereof must not be exceeded (see *"Chapter 11, Technical Data, on page 105"*).
- Do not attach or mount any additional loads to the extension arm, spring arms, adaptations and end devices.

Collision damage

In case of collision with other appliances, walls or ceilings, the pendant system can be damaged and important patient care systems can fail:

- After a collision, the pendant system must be inspected for damage.
- In case of doubt, contact your supplier.

2.5.2 Mounting / dismantling

NOTICE**Electric shock hazard**

To prevent the risk of electric shock, the appliance may only be connected to a supply network equipped with a protective conductor.

- The appliance must be connected in such a way that it can be disconnected from the mains at all poles and at the same time.

WARNING**Electric shock hazard**

Power supply cables may have been laid in the pendant system. Contact with energised components presents a danger to life from electric shock. Prior to mounting/dismantling work, disconnect the appliance from the mains:

- Disconnect all the poles from the mains and prevent the appliance from being switched back on again.
- Check whether the end device (e.g. flat screen, OR lamp, etc.) is de-energised.

**Sudden release of spring arm**

When dismantling the end device from the spring arm [4] - [6], the spring arm [4] - [6] can suddenly jump up and may cause serious injury:

- Before removing the end device, make sure that the vertical lift of the spring arm is restricted to the lowermost end stop (horizontal position).

**Risk of parts falling off**

During all dismantling and mounting work, it must be ensured that nobody is located in the area underneath the pendant system.

2.5.3 Cleaning and disinfection

Cleaning

WARNING**Risk of contamination and infection of the patient**

Parts of the pendant system and the adaptations are made of plastic; solvents can dissolve plastic materials. Strong acids, bases and agents with an alcoholic strength of more than 60% can lead to the plastic materials becoming brittle. Damaged parts can fall into open wounds. Liquid cleaning agents penetrating the pendant system and the adaptations can drip into open wounds.

Disinfection

CAUTION**Health hazard**

Disinfectants can contain substances hazardous to health which, when in contact with the skin and eyes, can cause injuries or affect the respiratory organs when inhaled. Observe the protective measures:

- Observe the hygiene regulations.
- Adhere to the disinfectant manufacturer's instructions.
- Perform surface disinfection every working day and in case of contamination.

2.5.4 Maintenance work

WARNING

Electric shock hazard

Power supply cables are laid in the pendant system and the adaptations **7** - **9**. Contact with energised components presents a danger to life from electric shock. Disconnect the pendant system from the mains before any maintenance work:

- Disconnect all the poles of the pendant system from the mains and prevent it from being switched back on again.
- Check whether the end device (e.g. flat screen, OR lamp, etc.) is de-energised.
- Wait until the end device (e.g. flat screen, OR lamp, etc.) has cooled down.

CAUTION

Maintenance work

- The pendant system and the adaption **7** - **9** must be checked as described in *“Chapter 13, Inspection Plan, on page 108”*.
- In case of failure or damage, please contact your supplier.

Warranty

Ondal Medical Systems GmbH warrants the functional reliability of the pendant system only under the condition that:

- No structural alterations are made to the pendant system. Unauthorised modifications or conversions of the pendant system are not permitted for safety reasons.
- Only genuine spare parts or accessories and those defined and approved by Ondal Medical Systems are used. The use of other parts may involve unknown risks and must be avoided in all cases.
- Inspections and maintenance work must be carried out at the specified time intervals.

3.1 Parts and control elements

Please open the hinged side.

- 1 Canopy (ceiling-mounted version) / wall bearing (wall-mounted version)
- 2 Ceiling tube (ceiling-mounted version only)
Variable length for compensating different ceiling heights
- 3 Extension arm in different lengths
- 4 ACROBAT 2000 spring arm (standard version)
- 5 LCH spring arm (Low Ceiling Height)
- 6 NRH spring arm (Low Ceiling Height)

Approved accessories for the spring arm 4 (optional):
(see also "Chapter 12" on page 107)

- 7 Connecting piece – for mounting an end device (e.g. heat radiator)
- 8 Bracket – for mounting an end device (e.g. OR lamp, camera)
- 9 Monitor connecting piece – for mounting a flat screen
- 10 Vario View flat screen support, single connection – for mounting a flat screen
- 11 Vario View flat screen support, dual connection – for mounting 2 flat screens

Approved accessories for the NRH spring arm 5 / LCH spring arm 6 (optional):

- 8 Bracket

Only valid when combined with the Installation and Operating Instructions of the Vario View flat screen supports

3.2 Validity of these Operating Instructions

The pendant system and the **Vario View flat screen support** [10], [11] are described in separate Installation and Operating Instructions.

These Operating Instructions of the **pendant system** are only valid and complete when combined with the Installation and Operating Instructions of the **Vario View flat screen support** [10], [11]s.

For information on installation and operation of the **Vario View flat screen support** [10], [11] refer to the following documents:

- Installation Instructions: No. 1522529
- Operating Instructions: No. 1522530

3.3 Appliance description


Canopy / wall bearing	The points for connecting the electrical cables of the pendant system are located underneath the canopy / wall bearing [1]. The canopy / wall bearing [1] may only be removed by specialist personnel authorised by the operator.
Ceiling tube (ceiling-mounted version only)	The ceiling tube [2] (ceiling-mounted version only) compensates different ceiling heights in order to ensure that the adaption [7] - [9] with the end device (e.g. flat screen, OR lamp, etc.) is positioned at the working height desired.
Extension arm	The extension arms [3] with the internal supply cables (e.g. for a flat screen) can be swiveled 330 degrees, and with internal plug couplings e.g. for OR lamps >360 degrees.
Twist protection	In order to prevent the internal supply cables tearing off, the swivel range of the extension arm [3] is restricted.
Spring arm with internal supply cables	The spring arm [4] with the internal supply cables, e.g. for a flat screen, can be swiveled 330 degrees horizontally.
Twist protection	In order to prevent the internal supply cables tearing off, the swivel range of the spring arm [4] is restricted.
Spring arm with internal plug couplings	The spring arms [4] - [6] with internal plug couplings, e.g. for OR lamps, can be rotated >360 degrees.
	In addition, the spring arms [4] - [6] are height adjustable.
Vertical lift of the spring arms	In order to prevent the spring arm [4] - [6] hitting the ceiling, its vertical lift can be restricted accordingly. The vertical lift was defined during installation.
Adaptations	The adaptions [7] - [9] serve as a support for holding the end devices (e.g. flat screen, OR lamp, etc.) as specified in " <i>Chapter 12, Optional Accessories, on page 107</i> ".

3.4 Functional description

Extension arm	The extension arms [3] serves for the horizontal positioning of the adaption [7] - [9] with the end device (e.g. flat screen, OR lamp, etc.). The swivel motions can be restricted through internal limit stops.
Brakes on the extension arm	The extension arms [3] can be (optionally) equipped with friction brakes in order to ensure that the extension arms [3], the spring arm [4] - [6], the adaption [7] - [9] and the end device (e.g. flat screen, OR lamp, etc.) remain stable in any set position.
Brakes	The extension arms [3] and the spring arm [4] - [6] are equipped with brakes in order to ensure that they remain stable in any set position. For more detailed information on how to adjust the brake force, refer to <i>"Chapter 9.3" on page 97.</i>
Spring arm	The spring arm [4] - [6] serves for the horizontal and vertical positioning of the adaption [7] - [9] with the end device (e.g. flat screen, OR lamp, etc.).
Functioning of the springs in the spring arm	To facilitate the positioning of the adaption [7] - [9] with the end device (e.g. flat screen, OR lamp, etc.), a spring is installed in the spring arm [4] - [6]. The spring compensates the weight of the adaption [7] - [9] with the end device (e.g. flat screen, OR lamp, etc.) and is available in different spring equipment versions (see <i>"Chapter 8.2.3" on page 82</i> , <i>"Chapter 8.2.4" on page 83</i> and <i>"Chapter 8.2.5" on page 84</i>).
Spring tension adjustment on the spring arm	If the spring arm [4] - [6] with the end device (e.g. flat screen, OR lamp, etc.) moves down, or if a new end device has been mounted, the spring tension must be readjusted. For more detailed information on how to adjust the spring tension, refer to <i>"Chapter 9.2" on page 95.</i>
Vertical lift of the spring arms	In order to prevent collisions with the ceiling or other components, the vertical lift of the spring arms [4] - [6] can be restricted. The vertical lift was defined during installation. For more detailed information on how to adjust the vertical lift, refer to <i>"Chapter 9.4" on page 100.</i>
Adaptations and end devices	The adaption [7] - [9] have been designed for the transportation of the end device (e.g. flat screen, OR lamp, etc.).

Initial commissioning	<ol style="list-style-type: none">1. The pendant system must be installed. Installation instructions are included in the scope of delivery of the product.2. For commissioning following installation, proper initial commissioning must be carried out for the entire pendant system.
Functional test	<p>Prior to using the pendant system on a patient for the first time, a functional test at the place of installation must be carried out by the operator or a person authorised by the operator; the persons authorised by the operator must have been duly instructed by the operator.</p> <p>To use the pendant system for medical operation, these Operating Instructions only apply if:</p> <ol style="list-style-type: none">1. The functional reliability of the pendant system has been ensured.2. The maximum permissible loading capacity (payload) has been safely determined and is indicated on a label attached to the spring arm or adaptation.3. Proper functioning has been approved by the operator during initial commissioning and documented by signing a declaration of acceptance. <p>The following points must be observed during handover to the operator:</p> <ol style="list-style-type: none">1. The pendant system must be subjected to extensive tests before being handed over to the operator.2. Handover must be documented in writing including confirmation by the operator.3. On handover, the operator must be instructed in the functioning and effect of the maximum loading capacity (payload).4. In addition, the operator must be instructed in the functioning, operation, cleaning and disinfection of the pendant system during the handover procedure.5. Permitted adjustments in accordance with the Operating Instructions included in the scope of delivery must be explained to the operator on handover.6. On completion of the instruction, an instruction certificate must be created and signed in order to document that the operator/user has understood the special operator control actions required for normal use.

Figure 3: Positioning the pendant system

5.1 Positioning the pendant system

(See "Figure 3")

The **extension arms** [3], the **spring arms** [4] - [6]s and the **adaption** [7] - [9] with the end device (e.g. flat screen, OR lamp, etc.) can be positioned easily.

The swivel range and the vertical lift can be restricted through internal end stops.

5.1.1 Swiveling the pendant system

(See "Figure 3")

NOTICE

Damage to the pendant system

To prevent damage to the pendant system:

- Do not hit the end stops hard;
 - avoid collisions with other components or systems.
-
- Slowly swivel the end device (e.g. flat screen, OR lamp, etc.).
 - ➔ The swivel range ends at the internal end stops of the **extension arms** [3] and the **spring arms** [4] - [6].

5.1.2 Adjusting the height of the pendant system

NOTICE

Damage to the pendant system

If the end device (e.g. flat screen, OR lamp, etc.) is moved upwards, there is a risk that it collides with other components:

- Prior to adjusting the height, check for potential risks of collision.
-
- Slowly adjust the height of the end device (e.g. flat screen, OR lamp, etc.).
 - ➔ The height adjustment ends at the internal end stops of the **spring arms** [4] - [6]s.

If the components of the pendant system, the **spring arms** [4] - [6], the **adaption** [7] - [9] or the end device (e.g. flat screen, OR lamp, etc.) do not remain stable at the set position, the spring tension must be adjusted as described in "Chapter 9.2" on page 95 or the brake force must be adjusted as described in "Chapter 9.3" on page 97.

6.1 General safety instructions

WARNING

Electric shock hazard

The appliances can carry an electric current and must be treated with the utmost care during cleaning and disinfection.

- If a mains plug exists, pull the mains plug.
- Do not apply spray cleaning and/or spray disinfection.
- Do not spray liquid into power sockets, gas sockets or appliance openings and prevent the penetration of liquids.

6.2 Cleaning

Follow the safety instructions

1. Follow the general safety instructions prescribed in *"Chapter 6.1"*.

Cleaning agents

Recommended cleaning agents

Use a mild soap solution or a regular dishwashing product.

2. Wipe the surfaces of the appliances with a moderately moist cloth; add a mild soap solution (dishwashing product) if required.
3. Afterwards, carefully wipe the surfaces dry with a clean cloth.

6.3 Disinfection

Follow the safety instructions

1. Follow the general safety instructions prescribed in *"Chapter 6.1"*.

Disinfection method

Wiping disinfection is the standardised disinfection method prescribed for the pendant system. Hygiene regulations and related safety instructions for the disinfection methods to be applied must be defined by the operator.

- In case of contamination with potentially infectious material (e.g. blood, body secretion or excrement) the surfaces must be immediately and specifically disinfected.
- Make sure that you apply the disinfectant in the correct concentration.
- For surface disinfection do not spray, but wipe, the surfaces.
- Wiped surfaced may only be used after the disinfectant has dried.

WARNING

The appliance is not suitable for sterilisation

Avoid damage

Make sure that no liquid penetrates the appliance whilst cleaning it.

To prevent damage to plastic parts, refrain from using abrasives or alkaline, acidic or corrosive cleaning agents.

Do not use bleaching agents on stainless steel parts.

Use trained technical specialists only and abide by national regulations.

Cleaning/disinfection must be carried out by trained technical specialists only. The requirements of the national hygiene and disinfection committee must be complied with.

7.1 General safety instructions

⚠ WARNING

Electric shock hazard

Power supply cables are laid in the appliance. Contact with energised components presents a danger to life from electric shock. Disconnect the pendant system from the mains before any maintenance work:

- Disconnect all the poles of the pendant system from the mains and prevent it from being switched back on again.
- Check whether the end device (e.g. flat screen, OR lamp, etc.) is de-energised.
- Wait until the end device (e.g. flat screen, OR lamp, etc.) has cooled down.

Necessary maintenance work must be carried out in accordance with the inspection plan in "Chapter 13, Inspection Plan, on page 108".

7.2 Built-in components from third-party manufacturers

Built-in components from third-party manufacturers must be inspected and maintained as prescribed in the applicable Operating Instructions.

Third-party manufacturers

7.3 Maintenance of the securing segment


(See "Figure 4")

Checking the securing segment on the spring arms

Check that the securing segment of the **spring arms** [4] - [6] is at least 1.5 mm thick and 9.5 mm wide, and grease it; it must be replaced if the minimum thickness or width are not complied with.

- Dismantle and check the **securing segment** ① of the **spring arm** [4] as described in "Chapter 8.3" on page 85.
- Dismantle and check the **securing segment** ① of the **LCH spring arm** [6] as described in "Chapter 8.7" on page 91.

Figure 4: Maintenance of the securing

segment

7.4 Repeated inspection

Repeated inspections must be carried out in accordance with DIN EN 62353.

8.1 General safety instructions

WARNING

Electric shock hazard

Power supply cables are laid in the appliance. Contact with energised components presents a danger to life from electric shock. Prior to mounting and dismantling work, disconnect the appliance from the mains.

- Disconnect all the poles of the pendant system from the mains and prevent it from being switched back on again.
- Check whether the end device (e.g. flat screen, OR lamp, etc.) is de-energised.
- Wait until the end device (e.g. flat screen, OR lamp, etc.) has cooled down.

WARNING

Risk of parts falling off

During all mounting / dismantling work, it must be ensured that nobody is located in the area underneath the pendant system.

Pendant system crashing

The spring arm **4** - **6** and the extension arms **3** are adapted to each other with regard to their maximum loading capacity. The spring arm **4** - **6** must therefore not be dismantled from the extension arms **3** and mounted to another extension arms **3**.

CAUTION

Carrying out mounting and dismantling work

Mounting and dismantling work may only be carried out by specialist personnel authorised and instructed by the operator.

8.2 Calculating the maximum loading capacity (payload)

⚠ WARNING


Pendant system crashing

If the maximum loading capacity has been exceeded, there is a risk that the pendant system or components thereof may disengage from the fastening device and crash:

- The maximum loading capacity on the **spring arm** [4] - [6] must not be exceeded!
- On the **spring arm** [4] (standard version), the **adaption** [7] - [9] mounted and the end device (e.g. flat screen) may only introduce vertical weight forces, and no bending moments, into the **spring arm** [4] (standard version).

Figure 5: What is the maximum loading capacity (payload)?

8.2.1 What is the maximum loading capacity (payload)?

(See "Figure 5")

The maximum loading capacity (payload of e.g. flat screens, OR lamps) corresponds to the maximum weight which may be carried by a **spring arm** [4] - [6].

The various **spring arms** [4] - [6] can carry different loading capacities:

- For more detailed information on the maximum loading capacities of the **spring arms** [4] - [6], refer to "Chapter 11, Technical Data, on page 105".

Calculate the maximum loading capacity (payload) of the different **spring arms** [4] - [6]:



- **spring arm** [4] as described in "Chapter 8.2.3" on page 82.
- **NRH spring arm** [5] as described in "Chapter 8.2.4" on page 83.
- **LCH spring arm** [6] as described in "Chapter 8.2.5" on page 84.

8.2.2 Determining the maximum loading capacity (payload)

The total weight of the **adaption** [7] - [9] and the end device (e.g. flat screen, OR lamp, etc.) must not exceed the maximum loading capacity (payload) of the **spring arms** [4] - [6] and must therefore be calculated.

The maximum loading capacity (payload) of the **spring arms** [4] - [6] is indicated on the rating plate on the top side of the **spring arms** [4] - [6].

Figure 6: Calculating the maximum loading capacity of the ACROBAT 2000 spring arm

8.2.3 Calculating the maximum loading capacity (payload) of the ACROBAT 2000 spring arm

(See "Figure 6")

⚠ WARNING

Sudden release of spring arm

The spring arm (4) contains a powerful spring. When the end device is dismantled, the spring arm can suddenly jump up and may cause serious injury:

- Before removing the end device (e.g. OR lamp), make sure that the vertical lift of the spring arm (4) as described in "Chapter 9.4.1" on page 100 is restricted to the lowermost end stop (horizontal position).

Maximum loading capacity of the spring arm (standard version)

The maximum loading capacity indicated corresponds to the maximum weight which can be carried by a spring arm. For example, a maximum loading capacity of 21 kg* is indicated on the rating plate. To check whether the spring arm can bear the actual load of an end device, the actual load must first be calculated.

The total weight of the adaptation and the end device must not exceed the maximum loading capacity (spring equipment) of the spring arms.

Possible loading capacities (spring equipment versions)

- 1.5 - 3.5 kg; 3.5 - 7 kg; 7 - 12 kg; 9 - 15 kg; 12 - 18 kg; 13.5 - 21 kg.
- The spring equipment mounted is indicated on the rating plate on the top side of the spring arm (4).

Calculating the actual load

The dead weight of the adaption (7) - (9) mounted and the end device (e.g. OR lamp) must be subtracted from the maximum loading capacity of 21kg, and the result must be greater than 0.

* for spring equipment with a loading capacity of up to 21kg. For other spring equipment versions, the corresponding value must be used.

Configuration example

- Max. loading capacity of the spring arm (4) with spring equipment mounted = up to 21 kg*,
- dead weight of the VESA adapter (9) = 2.0kg

Calculation of the maximum loading capacity

- Max. loading capacity (payload) = 21kg* - 2.0kg = 19.0kg.

Evaluation

The weight of the end device to be mounted must not exceed 19.0 kg.

Maximum torque

The torque produced on the front joint (1) through the off-centre attachment of a load F and/or an actuating force F must not exceed 25 Nm.

Calculating the torque


Torque produced through manual force:

This torque is determined from the manual force F (N) and the vertical distance Lx (m).

Torque produced by the end device:

This torque is determined from the force of the weight of the end device F (N) and the adapter (N) and the horizontal distance Ly (m).

Figure 7: Calculating the maximum loading capacity (payload) of the NRH spring arm

8.2.4 Calculating the maximum loading capacity (payload) of the NRH spring arm

(See "Figure 7")

WARNING

Sudden release of spring arm

The NRH spring arm (5) contains a powerful spring. When the end device is dismantled, the spring arm can suddenly jump up and may cause serious injury:

- Before removing the end device (e.g. OR lamp), make sure that the vertical lift of the NRH spring arm (5) as described in "Chapter 9.4.2" on page 100 is restricted to the lowermost end stop (horizontal position).

The weight of an end device (e.g. OR lamp) acts as a lever on the spring arm from the spring arm axis (1) up to the centre of gravity of the end device (2).

Maximum loading capacity of the NRH spring arm

The maximum loading capacity indicated corresponds to the maximum weight which can be borne by a spring arm. For example, a maximum loading capacity of 176 Nm is indicated on the rating plate. To check whether the spring arm can bear the actual load of an end device, the actual load must first be calculated.

The total weight of the adaptation and the end device must not exceed the maximum loading capacity (spring equipment) of the spring arms.

Possible loading capacities (spring equipment versions)

- 95 - 160Nm; 95 - 176 Nm
- The spring equipment mounted is indicated on the rating plate on the top side of the NRH spring arm (5).

Calculating the actual load

The actual load is calculated from the distance from the joint (1) of the spring arm to the centre of gravity of the end device (2) multiplied by the weight of the adaptation and the end device in Nm (1 kg = 9.81 N).

- The maximum loading capacity (payload) of the NRH spring arm (5) is 176Nm*.
- * for spring equipment with a loading capacity of up to 176Nm. For other spring equipment versions, the corresponding value must be used.

Configuration example

- Maximum loading capacity of the NRH spring arm (5) = 95 - 176 Nm*
- Length of the NRH spring arm (5) = 745 mm = 0.745m
- Distance l between the centre of gravity of the appliance and the NRH spring arm (5) = 300mm = 0.3m
- Dead weight of the adaptation (8) (bracket for OR lamp) = 5kg
- Dead weight of the end device (e.g. OR lamp) = 12.0 kg


Calculation of the maximum loading capacity

- Total length of the lever = 0.745 m + 0.3 m = 1.045m
- Adaptation weight (bracket) in N = 5kg x 9.81N = 49.05N
- Weight of the end device (e.g. OR lamp) in N = 12kg x 9.81 N = 117.72N
- Maximum loading capacity = 1.045m x (49.05N + 117.72N) = 174.28 Nm

Evaluation

The calculated weight of the OR lamp with bracket (8) to be mounted of 174.28 Nm is below the maximum load bearing capacity of the NRH spring arm (5) of 176 Nm.

Figure 8: Calculating the maximum loading capacity (payload) of the LCH spring arm

8.2.5 Calculating the maximum loading capacity (payload) of the LCH spring arm

(See "Figure 8")

⚠ WARNING

Sudden release of spring arm

The LCH spring arm (6) contains a powerful spring. When the end device is dismantled, the spring arm can suddenly jump up and may cause serious injury:

- Before removing the end device (e.g. OR lamp), make sure that the vertical lift of the spring arm (4) as described in "Chapter 9.4.3" on page 101 is restricted to the lowermost end stop (horizontal position).

The weight of an end device (e.g. OR lamp) acts as a lever on the spring arm from the spring arm joint (1) up to the centre of gravity of the end device (2).

Maximum loading capacity of the LCH spring arm

The maximum loading capacity indicated corresponds to the maximum weight which can be borne by a spring arm. For example, a maximum loading capacity of 176 Nm is indicated on the rating plate. To check whether the spring arm can bear the actual load of an end device, the actual load must first be calculated.

The total weight of the adaptation and the end device must not exceed the maximum loading capacity (spring equipment) of the spring arms.

Possible loading capacities (spring equipment versions)

- 95 - 160Nm; 95 - 176Nm
- The spring equipment mounted is indicated on the rating plate on the top side of the LCH spring arm (6).

Calculating the actual load

The actual load is calculated from the distance from the joint (1) of the spring arm to the centre of gravity of the end device (2) multiplied by the weight of the adaptation and the end device in Nm (1 kg = 9.81N).

- The maximum loading capacity (payload) of the LCH spring arm (6) is 176Nm*.
- * for spring equipment with a loading capacity of up to 176Nm. For other spring equipment versions, the corresponding value must be used.

Configuration example

- Maximum load bearing capacity of the spring arm LCH = 176 Nm
- Length of the spring arm = 700 mm = 0.7 m
- Distance I between the centre of gravity of the appliance and the spring arm = 300 mm = 0.3 m
- dead weight of the adaptation (8) and the end device (e.g. OR lamp) = 16 kg


Calculation of the maximum loading capacity

- Total length of the lever = 0.7 m + 0.3 m = 1.0 m
- Adaptation weight (bracket) in Nm = 5kg x 9.81N = 49.05N
- Weight of the end device (e.g. OR lamp) in Nm = 12kg x 9.81 N = 117.72N
- Maximum loading capacity = 1.0m x (49.05N + 117.72N) = 166.77Nm

Evaluation

The calculated weight of the OR lamp with bracket (8) to be mounted of 166.77 Nm is below the maximum load bearing capacity of the LCH spring arm (6) of 176 Nm.

Figure 9: Dismantling the end device from the spring arm (standard version)

8.3 Dismantling the end device from the ACROBAT2000 spring arm

(See "Figure 9")

Follow the safety instructions

Follow the general safety instructions prescribed in "Chapter 8.1" on page 80.

1. Check whether the end device (e.g. OR lamp) on the **spring arm** [4] can now be rotated more than 360 degrees.
- If this is the case, an **internal 3-, 5- or 7/9-pole plug coupling** [4]-[6] has been mounted and the end device (e.g. OR lamp) may be dismantled by you as a technician authorised by the operator.

⚠ WARNING

Sudden release of spring arm

The **spring arm** [4] contains a powerful spring. When the end device (e.g. OR lamp) is dismantled, the spring arm can suddenly jump up and may cause serious injury:

- Before removing the end device (e.g. OR lamp), make sure that the vertical lift of the **spring arm** [4] as described in "Chapter 9.4.1" on page 100 is restricted to the lowermost end stop (horizontal position).
2. Restrict the vertical lift of the **spring arm** [4] to the horizontal position as described in "Chapter 9.4.1" on page 100.
 - ➔ It must no longer be possible to move the **spring arm** [4] upwards beyond the horizontal position.
 3. Unscrew the **sleeve securing screw M3 x 8mm** [2] and push the **sleeve** [1] upwards.
 4. To relieve the **securing segment** [3] slightly lift the end device and remove it using a suitable slotted screwdriver.

NOTE – Different plug couplings

Depending on the version of the end device (e.g. OR lamp), the **spring arm** [4] is equipped with a 3-, 5- or 7/9-pole plug coupling.

For the 3-pole [4] and 5-pole plug couplings [5] only:

- Carefully remove the end device (e.g. OR lamp) from the **spring arm** [4], keeping it straightly aligned.

For the 7/9-pole plug coupling [6] only:

- Carefully remove the end device (e.g. OR lamp) from the **spring arm** [4], keeping it straightly aligned.
- Screw the 3 **plug adapters** [8] apart, disengage the **6-pole plug** [7] and remove the end device (e.g. OR lamp).

5. Insert the **securing segment** [3], push the **sleeve** [1] downwards and screw the **sleeve securing screw M3 x 8mm** [2] back in.

Marking

If you do not mount another end device (e.g. OR lamp) directly after dismantling, the plug coupling in the **spring arm** [4] is energised whilst the mains supply is switched on.

1. Mark the mains supply in the room with a warning sign:
ELECTRIC SHOCK HAZARD – DO NOT SWITCH ON THE POWER SUPPLY!
2. Mount the end device (e.g. OR lamp) later as described in "Chapter 8.4" on page 86.

8.4 Mounting the end device to the ACROBAT2000 spring arm

Follow the safety instructions

Follow the general safety instructions prescribed in "Chapter 8.1" on page 80.

Maximum loading capacity

Calculate the maximum loading capacity as described in "Chapter 8.2.3" on page 82.

Preparing the installation

(See "Figure 10")

1. Unscrew the **sleeve securing screw M3 x 8mm** ② and push the **sleeve** ① upwards.
2. Screw the **sleeve securing screw M3 x 8mm** ② back in in order to fasten the **sleeve** ①.
3. Remove the **securing segment** ⑤ using a suitable slotted screwdriver.

Only for version with 3-pole plug coupling

(See "Figure 10")

1. Grease the pivot groove and the **pivot** ④ using a suitable friction bearing grease (e.g. Microgleit GP 360).
2. Position the **pivot with 3-pole plug coupling** ④, keeping it straightly aligned, and establish the **assembly connection** ③ by gently pushing towards the **spring arm** ④.

Only for version with 5-pole plug coupling

(See "Figure 11")

1. Grease the pivot groove and the **pivot** ④ using a suitable friction bearing grease (e.g. Microgleit GP 360).

NOTE – Plug positioning through the grooves

- The **5-pole plug coupling** ④ has a **narrow groove** ⑥ and a **wide groove** ⑤.
 - To establish the **plug connection** ③, the **grooves** ⑤/⑥ must fit into the corresponding grooves of the counter piece.
2. Position the **pivot with 5-pole plug coupling** ④, keeping it straightly aligned, and establish the **assembly connection** ③ by gently pushing towards the **spring arm** ④.

Figure 10: Version with 3-pole plug coupling


Figure 11: Version with 5-pole plug coupling



Figure 12: Version with 7/9-pole plug coupling

Only for version with 7/9-pole plug coupling

(See "Figure 12")

Connecting the plugs

1. Connect the **6-pole plug** (4).
2. Screw 3 **plug adapters** (5) into each other in accordance with the colour code.
3. Gently push the **6-pole plug** (4) and the 3 **plug adapters** (5) into the **spring arm** (4).

Mounting the end device

4. Grease the pivot groove and the **pivot** (3) using a suitable friction bearing grease (e.g. Microgleit GP 360).
5. Position the **pivot** (3), keeping it straightly aligned, and establish the **assembly connection** (3) by gently pushing towards the **spring arm** (4).

Completing the installation

(See "Figure 12")

⚠ WARNING

End device crashing

When using an incorrect securing segment (5) the end device (e.g. OR lamp) can crash and the spring arm (4) can suddenly jump upwards:

- Only use a **genuine securing segment** (5).

1. Insert the **genuine securing segment** (5), unscrew the **sleeve securing screw M3 x 8mm** (2) and push the **sleeve** (1) downwards.

⚠ WARNING

End device crashing

The sleeve (1) secures the securing segment (5) and thus ensures that the end device in the spring arm (4) is securely in place:

- Check that the **sleeve** (1) covers the **securing segment** (3).
- Secure the **sleeve** (1) using the **sleeve securing screw M3 x 8mm** (2).

2. Tighten the **sleeve securing screw M3 x 8mm** (2) again.
3. Check that the end device (e.g. OR lamp) is securely in place.

Extending the vertical lift

4. If required, extend the vertical lift of the **spring arm** (4) as described in "*Chapter 9.4.1*" on page 100.

Adjustments

5. Adjust the brake force and spring tension as described in "*Chapter 9*" on page 95.

Function test


6. Connect the system to the mains.
7. Perform a function test.

ACROBAT 2000 appliances are only intended for the connection of luminaires for medical diagnosis, surgical luminaires and flat screens bearing the CE mark and with a maximum payload of 21 kg.

For technical assistance for creating a medical-electrical (ME) system refer to DIN EN 60601-1, Part 16.

The manufacturer / marketer of the ME system must provide a declaration in accordance with Article 12 of 93/42/EEC (Medical Device Directive, MDD).

Figure 13: Dismantling the end device from the NRH spring arm

8.5 Dismantling the end device from the NRH spring arm

(See "Figure 13")

Follow the safety instructions

Follow the general safety instructions prescribed in "Chapter 8.1" on page 80.

1. Check whether the end device (e.g. OR lamp) on the NRH spring arm [5] can be rotated more than 360 degrees.
2. If this is the case, an **internal 3- or 5-pole plug coupling** [5], [7] has been mounted and the end device (e.g. OR lamp) may be dismantled by you as a technician authorised by the operator.

⚠ WARNING

Sudden release of spring arm

The NRH spring arm [5] contains a powerful spring. When the end device (e.g. OR lamp) is dismantled, the spring arm can suddenly jump up and may cause serious injury:

- Before removing the end device (e.g. OR lamp), make sure that the vertical lift of the NRH spring arm [5] as described in "Chapter 9.4.2" on page 100 is restricted to the lowermost end stop (horizontal position).

3. Restrict the NRH spring arm [5] to the horizontal position as described in "Chapter 9.4.2" on page 100.
- ➔ It must no longer be possible to move the NRH spring arm [5] upwards beyond the horizontal position.

Dismantling the brake and sleeve securing screws

4. Turn the fitting aperture [3] of the sleeve [1] downwards.
5. Unscrew the brake screw [6] at the bottom of the sleeve [1].
6. Rotate the sleeve [1] 90 degrees and unscrew the first end device securing screw [2].
7. Rotate the sleeve [1] 180 degrees and unscrew the second end device securing screw [2].

NOTE – Different plug couplings

Depending on the version of the end device (e.g. OR lamp), the NRH spring arm [5] is equipped with a **3-pole** [5] or **5-pole plug coupling** [7].


8. Carefully remove the end device (e.g. OR lamp) from the NRH spring arm [5], keeping it precisely axially aligned.
9. Turn the fitting aperture [3] of the sleeve [1] downwards.
10. Screw in the brake screw [6] at the bottom of the sleeve [1].
11. Rotate the sleeve [1] 90 degrees and screw in the first end device securing screw [2].
12. Rotate the sleeve [1] 180 degrees and screw in the second end device securing screw [2].

Marking

If you do not mount another end device (e.g. OR lamp) directly after dismantling, the plug coupling in the NRH spring arm [5] is energised whilst the mains supply is switched on.

1. Mark the mains supply in the room with a warning sign:
ELECTRIC SHOCK HAZARD – DO NOT SWITCH ON THE POWER SUPPLY!
2. Mount the end device (e.g. OR lamp) later as described in "Chapter 8.6" on page 89.

Figure 14: Version with 3-pole plug coupling

8.6 Mounting the end device to the NRH spring arm

Follow the safety instructions

Follow the general safety instructions prescribed in "Chapter 8.1" on page 80.

Maximum loading capacity

Calculate the maximum loading capacity (payload) as described in "Chapter 8.2" on page 81.

Preparing the installation

(See "Figure 14")

1. Turn the **fitting aperture** ③ of the **sleeve** ① downwards.
2. Unscrew the **brake screw** ⑥ at the bottom of the **sleeve** ①.
3. Rotate the **sleeve** ① 90 degrees and unscrew the first **end device securing screw** ②.
4. Rotate the **sleeve** ① 180 degrees and unscrew the second **end device securing screw** ②.

Only for version with 3-pole plug coupling

(See "Figure 14")


5. Grease the pivot groove and the **pivot** ⑤ using a suitable anti-friction bearing grease (e.g. Optimol).
6. Position the **pivot with 3-pole plug coupling** ⑤, keeping it straightly aligned, and establish the **assembly connection** ④ by gently pushing towards the **NRH spring arm** ⑤.

ACROBAT 2000 appliances are only intended for the connection of luminaires for medical diagnosis, surgical luminaires and flat screens bearing the CE mark and with a maximum payload of 21 kg.

For technical assistance for creating a medical-electrical (ME) system refer to DIN EN 60601-1, Part 16.

The manufacturer / marketer of the ME system must provide a declaration in accordance with Article 12 of 93/42/EEC (Medical Device Directive, MDD).

Figure 15: Version with 5-pole plug coupling

Only for version with 5-pole plug coupling

(See "Figure 15")

1. Grease the pivot groove and the pivot (5) using a suitable anti-friction bearing grease (e.g. Optimol).

NOTE – Plug positioning through the grooves

- The 5-pole plug coupling (5) has a narrow groove (7) and a wide groove (8).
 - To establish the plug connection (4), the grooves (7)/(8) must fit into the corresponding grooves of the counter piece.
2. Position the pivot with 5-pole plug coupling (5), keeping it straightly aligned, and establish the assembly connection (5) by gently pushing towards the NRH spring arm (5).

Completing the installation

(See "Figure 15")

1. Turn the fitting aperture (3) of the sleeve (1) downwards.
2. Screw in the brake screw (6) at the bottom of the sleeve (1).

⚠ WARNING

End device crashing

The sleeve (1) with the two end device securing screws (2) ensures that the end device is securely in place in the NRH spring arm (5):

- Check that the two end device securing screws (2) have been mounted and that the sleeve (1) covers the two end device securing screws (2).
3. Rotate the sleeve (1) 90 degrees and screw in the first end device securing screw (2).
 4. Rotate the sleeve (1) 180 degrees and screw in the second end device securing screw (2).
 5. Check that the end device (e.g. OR lamp) is securely in place.

Extending the vertical lift

6. If required, extend the vertical lift of the NRH spring arm (5) as described in "Chapter 9.4.2" on page 100.


Adjustments

7. Adjust the brake force and spring tension as described in "Chapter 9" on page 95.

Function test

8. Connect the system to the mains.
9. Perform a function test.

Figure 16: Dismantling the end device from the LCH spring arm

8.7 Dismantling the end device from the LCH spring arm

(See "Figure 16")

Follow the safety instructions

Follow the general safety instructions prescribed in "Chapter 8.1" on page 80.

1. Check whether the end device (e.g. OR lamp) on the LCH spring arm [6] can be rotated more than 360 degrees.
 - If this is the case, an **internal 3-, 5- or 7/9-pole plug coupling** [9]-[11] has been mounted and the end device (e.g. OR lamp) may be dismantled by you as a technician authorised by the operator.

Dismantling the cover panels

1. Unscrew 2 countersunk Phillips screws M3 x 6 mm [2].
2. Gently disengage each cover panel [1] on the two staggered straps [3] using a suitable slotted screwdriver.
3. Remove the 2 cover panels [1].

⚠ WARNING

Sudden release of spring arm

The LCH spring arm [6] contains a powerful spring. When the end device (e.g. OR lamp) is dismantled, the spring arm can suddenly jump up and may cause serious injury:

- Before removing the end device (e.g. OR lamp), make sure that the vertical lift of the LCH spring arm [6] as described in "Chapter 9.4.3" on page 101 is restricted to the lowermost end stop (horizontal position).

4. Restrict the LCH spring arm [6] to the horizontal position as described in "Chapter 9.4.3" on page 101.
 - ➔ It must no longer be possible to move the LCH spring arm [6] upwards beyond the horizontal position.
5. Slightly unscrew the 2 brake screws [7].
6. Unscrew the Allen countersunk screw M4 x 4 mm [4] and push the securing ring for the securing segment [6] backwards towards the spring arm.
7. To relieve the securing segment [5] slightly lift the end device and remove it using a suitable slotted screwdriver.

NOTE – Different plug couplings

Depending on the version of the end device (e.g. flat screen, OR lamp, etc.), the LCH spring arm [6] is equipped with a 3-, 5- or 7/9-pole plug coupling.

For the 3-pole [9] and 5-pole plug coupling [10] only:


- Disengage the assembly connection [8] and gently remove the adaptation with the end device (e.g. OR lamp) from the LCH spring arm [6], keeping it straightly axially aligned.

For the 7/9-pole plug coupling [11] only:

- Gently withdraw the adaptation with the end device (e.g. OR lamp) from the LCH spring arm [6].
- Screw the 3 plug adapters [13] apart, disengage the 6-pole plug [12] and remove the end device (e.g. OR lamp).

8. Insert the securing segment [5], push the securing ring for the securing segment [6] forwards and screw in the Allen countersunk screw M4 x 4 mm [4].

Figure 17: Mounting the cover panels

**Mounting the cover panels***(See "Figure 17")*

9. Position the **cover panel** ① with the internal **guide** ⑤ in such a way that the **guide** ⑤ protrudes into the **adjustment hole** ④.
10. Gently push the **cover panel** ① and fasten it to the spring arm using the **countersunk Phillips screw M3 x 6 mm** ②.
11. Gently push the second **cover panel** ① so that the 4 **straps** ③ snap into each other.
12. Screw the **cover panel** ① to the spring arm using the **countersunk Phillips screw M3 x 6 mm** ②.
13. Check that the **cover panels** ① are securely in place.
 - The **cover panels** ① must sit tightly on the spring arm without any gaps.

Marking

If you do not mount another end device (e.g. OR lamp) directly after dismantling, the plug coupling in the **LCH spring arm** ⑥ is energised whilst the mains supply is switched on.

1. Mark the mains supply in the room with a warning sign:
ELECTRIC SHOCK HAZARD – DO NOT SWITCH ON THE POWER SUPPLY!
2. Mount the end device (e.g. OR lamp) later as described in *"Chapter 8.8" on page 92*.

8.8 Mounting the end device to the LCH spring arm**Follow the safety instructions**

Follow the general safety instructions prescribed in *"Chapter 8.1" on page 80*.

Maximum loading capacity

Calculate the maximum loading capacity as described in *"Chapter 8.2" on page 81*.

Dismantling the cover panels

1. Dismantle the cover panels as described in *"Chapter 8.7" on page 91*.

Preparing the installation*(See "Figure 18")*

2. Unscrew the **Allen countersunk screw M4 x 4 mm** ① and push the **securing ring for the securing segment** ③ backwards.
3. Remove the **securing segment** ② using a suitable slotted screwdriver.
4. Slightly unscrew the 2 **brake screws** ④ so that they no longer protrude into the fitting aperture of the spring arm.

NOTE – Different plug couplings

Depending on the version of the end device (e.g. flat screen, OR lamp, etc.), the **LCH spring arm** ⑥ is equipped with a 3-, 5- or 7/9-pole plug coupling.

Only for version with 3-pole plug coupling*(See "Figure 18")***NOTICE****Damage to the adaptation / end device**

Grease the pivot and the pivot groove before mounting the end device.

1. Grease the pivot groove and the **pivot** ⑥ using a suitable anti-friction bearing grease (e.g. Optimol).

Figure 18: Version with 3-pole plug coupling


Figure 19: Version with 5-pole plug coupling



Figure 20: Version with 7/9-pole plug coupling

2. Position the **pivot with 3-pole plug coupling** ⑥, keeping it straightly aligned, and establish the **assembly connection** ⑤ by gently pushing towards the **LCH spring arm** ⑥.

Only for version with 5-pole plug coupling

(See "Figure 19")

1. Grease the pivot groove and the **pivot** ⑥ using a suitable anti-friction bearing grease (e.g. Optimol).

NOTE – Plug positioning through the grooves

- The **5-pole plug coupling** ⑥ has a **narrow groove** ⑦ and a **wide groove** ⑧.
- To establish the **plug connection** ⑤, the **grooves** ⑦/⑧ must fit into the corresponding grooves of the counter piece.

2. Position the **pivot with 5-pole plug coupling** ⑥, keeping it straightly aligned, and establish the **assembly connection** ⑤ by gently pushing towards the **LCH spring arm** ⑥.

Only for version with 7/9-pole plug coupling

(See "Figure 20")


Connecting the plug

1. Connect the **6-pole plug** ⑦.
2. Screw 3 **plug adapters** ⑧ into each other in accordance with the colour code.
3. Gently push the **6-pole plug** ⑦ and the plug adapters ⑧ into the **LCH spring arm** ⑥.

Mounting the end device

4. Grease the pivot groove and the **pivot** ⑥ using a suitable anti-friction bearing grease (e.g. Optimol).
5. Insert the **pivot** ⑥ into the **LCH spring arm** ⑥, keeping it precisely axially aligned.

Figure 21: Completing the installation

**Completing the installation***(See "Figure 21")***⚠ WARNING****End device crashing**

When inserting an incorrect securing segment ②, the adaptation / end device can crash and the LCH spring arm ⑥ can suddenly jump up:

- Only use a genuine securing segment ②.

The securing ring ③ secures the securing segment ② and thus ensures that the end device is securely in place in the LCH spring arm ⑥:

- The securing ring ③ must cover the securing segment ②.
- The Allen countersunk screw M4 x 4mm ① must be screwed in and tightened.

1. Insert a **genuine securing segment** ②.
2. Push the **securing ring for the securing segment** ③ forwards and then screw in and tighten the **Allen countersunk screw M4 x 4 mm** ① through the **securing ring for the securing segment** ③.
3. Check that the **pivot** ⑥ of the end device (e.g. flat screen, OR lamp, etc.) is securely in place.
 - The **pivot** ⑥ must sit tightly in the spring arm.
 - The **pivot** ⑥ must rotate easily without resistance.
4. Slightly tighten the opposite **brake screws** ④.

Extending the vertical lift

5. If required, extend the vertical lift of the **LCH spring arm** ⑥ as described in *"Chapter 9" on page 95*.

Adjustments

6. Adjust the brake force and spring tension as described in *"Chapter 9" on page 95*.

Mounting the cover panels

7. Mount the cover panels as described in *"Chapter 8.7" on page 91*.

Function test

8. Connect the system to the mains.
9. Perform a function test.

ACROBAT 2000 appliances are only intended for the connection of luminaires for medical diagnosis, surgical luminaires and flat screens bearing the CE mark and with a maximum payload of 21 kg.

For technical assistance for creating a medical-electrical (ME) system refer to DIN EN 60601-1, Part 16.

The manufacturer / marketer of the ME system must provide a declaration in accordance with Article 12 of 93/42/EEC (Medical Device Directive, MDD).

9.1 General safety instructions

⚠ WARNING

Electric shock hazard

Power supply cables are laid in the pendant system and the adaption [7] - [9]. Contact with energised components presents a danger to life from electric shock. Disconnect the pendant system from the mains before any adjustment work:

- Disconnect all the poles of the pendant system from the mains and prevent it from being switched back on again.
- Check whether the end device (e.g. flat screen, OR lamp, etc.) is de-energised.
- Wait until the end device (e.g. flat screen, OR lamp, etc.) has cooled down.

⚠ CAUTION

Performing adjustment work

Adjustment work may only be carried out by the operator.

9.2 Adjusting the spring tension

The **spring arm** [4] - [6] is equipped with a spring which compensates the weight of the adaptation and the end device.

Adjust the spring tension so that the **spring arm** [4] - [6] with the **adaption** [7] - [9] and the end device (e.g. flat screen, OR lamp, etc.) remains stable in any position.

If the **spring arm** [4] - [6] with the **adaption** [7] - [9] and the end device (e.g. flat screen, OR lamp, etc.) does not remain in position after adjusting the spring tension, the **spring arm** [4] - [6] must be replaced by a service technician.

Follow the safety instructions

Follow the general safety instructions prescribed in "Chapter 9.1" on page 95.

9.2.1 Adjusting the spring tension on the ACROBAT2000 and NRH spring arm (See "Figure 22")

Tool to be used

Use an Allen key (size 5).

Possible spring equipment on the ACROBAT 2000 spring arm

1.5 - 3.5kg; 3.5 - 7kg; 7 - 12kg; 9 - 15kg; 12 - 18kg; 13.5 - 21kg. The spring equipment mounted is indicated on the rating plate on the top side of the **spring arm** [4].

Possible spring equipment on the ACROBAT 2000 NRH spring arm

95 - 160Nm; 95 - 176Nm: The spring equipment mounted is indicated on the rating plate on the top side of the **NRH spring arm** [5].

1. Slightly pull down the **spring arm** [4], [5] in order to relieve the adjustment screw in the spring arm.
2. Insert the Allen key into the **adjustment hole** ①.

If the **spring arm** [4], [5] moves down, the spring tension is too low.

- Turn the Allen key to the left (anti-clockwise) as illustrated in the Figure.

If the **spring arm** [4], [5] moves up, the spring tension is too high.

- Turn the Allen key to the right (clockwise) as illustrated in the Figure.

3. Perform a function test.

Figure 22: Adjusting the spring tension on the ACROBAT 2000 and NRH spring arm




Figure 23: Adjusting the spring tension on the ACROBAT LCH spring arm

9.2.2 Adjusting the spring tension on the ACROBAT LCH spring arm

(See "Figure 23")

Tool to be used

Use an Allen key (size 5).

Possible spring equipment versions

95 - 160Nm; 95 - 176Nm: The spring equipment mounted is indicated on the rating plate on the top side of the **LCH spring arm** [6].

1. Slightly lift the **LCH spring arm** [6] in order to relieve the adjustment screw in the spring arm.
2. Insert the Allen key into the **adjustment hole** [1].

If the LCH spring arm [6] moves down, the spring tension is too low.

- Turn the Allen key to the left (anti-clockwise) as illustrated in the Figure.

If the LCH spring arm [6] moves up, the spring tension is too high.

- Turn the Allen key to the right (clockwise) as illustrated in the Figure.

3. Perform a function test.

9.3 Adjusting the brake force

The **extension arms** [3], the **NRH spring arm** [5] and the **LCH spring arm** [6] are equipped with optional brakes. The brakes hold the pendant system and the end device (e.g. flat screen, OR lamp, etc.) in the set position.

Adjust the brake force in such a way that the end device (e.g. flat screen, OR lamp, etc.) remains in the set position and can still be adjusted easily.


Follow the safety instructions

- Follow the general safety instructions in "*Chapter 9.1*" on page 95.

9.3.1 Adjusting the brake force on the pendant system

(See "Figure 24")

Figure 24: Adjusting the brake force on the pendant system

⚠ WARNING

Pendant system crashing

The **brake screws** [A], [B] are designed as **slotted screws**:

- No other screws must be loosened.

Tool to be used

Use a suitable slotted screwdriver.

To increase the brake force


- Insert the slotted screwdriver into the **brake screws** [A], [B] and alternately turn with the same number of revolutions in the clockwise (+) direction as illustrated in the Figure.

To reduce the brake force

- Insert the slotted screwdriver into the **brake screws** [A], [B] and alternately turn with the same number of revolutions in the anti-clockwise (-) direction.

Perform a function test.

Figure 25: Adjusting the brake force on the ACROBAT 2000 NRH spring arm

9.3.2 Adjusting the brake force on the ACROBAT2000 NRH spring arm
(See "Figure 25")

WARNING

Pendant system crashing

The brake screw ③ is designed as a slotted screw:

- No other screws must be loosened.

Tool to be used

Use a suitable slotted screwdriver.

1. Turn the fitting aperture ② in the cover ① downwards.
 - The brake screw ③ is located at the bottom of the connection.

To increase the brake force


- Insert the slotted screwdriver into the brake screw ③ and alternately turn with the same number of revolutions in the clockwise (+) direction as illustrated in the Figure ④.

To reduce the brake force

- Insert the slotted screwdriver into the brake screw ③ and alternately turn with the same number of revolutions in the anti-clockwise direction (-) as illustrated in the Figure ④.

2. Perform a function test.

Figure 26: Adjusting the brake force on the ACROBAT LCH spring arm

9.3.3 Adjusting the brake force on the ACROBAT LCH spring arm

(See "Figure 26")

Dismantling the cover panels

1. Unscrew the 2 countersunk Phillips screws M3 x 12 mm ②.
2. Gently disengage each cover panel ⑥ on the 2 staggered straps ⑤ using a suitable slotted screwdriver.
3. Remove the 2 cover panels ⑥.

Adjusting the brakes

⚠ WARNING

Pendant system crashing

The brake screws ③ are designed as slotted screws:

- All other screws on the spring arm must not be loosened.

Tool to be used

Use a suitable slotted screwdriver.

To increase the brake force

- Insert the slotted screwdriver into the brake screw ④ and alternately turn with the same number of revolutions in the clockwise (+) direction as illustrated in the Figure.

To reduce the brake force

- Insert the slotted screwdriver into the brake screw ④ and alternately turn with the same number of revolutions in the anti-clockwise direction (-) as illustrated in the Figure.

Mounting the cover panels

4. Position the cover ⑥ with the internal guide ③ in such a way that the guide ③ protrudes into the adjustment hole ①.
5. Gently push the cover panel ⑥ and fasten it to the spring arm using the countersunk Phillips screw M3 x 12mm ②.
6. Gently push the second cover ⑥ so that the 4 straps ⑤ snap into each other.
7. Screw the cover panel ⑥ to the spring arm tightly using the countersunk Phillips screw M3 x 12 mm ②.
8. Check that the cover panels ⑥ are securely in place.
 - The cover panels ⑥ must sit tightly on the spring arm without any gaps.
9. Perform a function test.

Figure 27: Adjusting the vertical lift on the ACROBAT 2000 spring arm





Figure 28: Adjusting the vertical lift on the ACROBAT 2000 NRH spring arm

9.4 Adjusting the vertical lift

The vertical lift prevents the **spring arms** 4 - 6 colliding with other devices or ceilings. Adjust the vertical lift so that collisions are excluded (the vertical lift can be restricted down to the horizontal position).

Follow the safety instructions

- Follow the general safety instructions in "Chapter 9.1" on page 95.

9.4.1 Adjusting the vertical lift on the ACROBAT 2000 spring arm

(See "Figure 27")

Tool to be used

Use an Allen key (size 5).

Adjusting the swiveling range

1. Insert the Allen key into the **adjustment hole** 1.
 - Slightly pull down the spring arm in order to relieve the adjustment screw in the spring arm.

To reduce the swivel range

- Turn the Allen key to the right (clockwise) as illustrated in the Figure.

To extend the swivel range

- Turn the Allen key to the left (anti-clockwise) as illustrated in the Figure.

2. Perform a function test.

9.4.2 Adjusting the vertical lift on the ACROBAT 2000 NRH spring arm

(See "Figure 28")

Tool to be used

Use the pin (Ø4mm x 110mm) available as an accessory included in the scope of delivery.

Dismantling the cover panels

1. Unscrew 2 **crosshead screws** 2 and remove the 2 **cover panels** 1.

Adjusting the vertical lift

2. Insert the pin (Ø 4mm x 110mm) into the perforated **adjusting nut** 4 in the **adjusting aperture** 3.

To reduce the swivel range

- Turn the **adjusting nut** 4 downwards in the anti-clockwise (-) direction as illustrated in the Figure.

To extend the swivel range



- Turn the **adjusting nut** 4 upwards in the clockwise (+) direction as illustrated in the Figure.

Mounting the cover panels

3. Position the **cover panels** 1 and screw in two **crosshead screws** 2.
4. Check that the **cover panels** 1 are securely in place.
 - The **cover panels** 1 must sit tightly on the spring arm without any gaps.

5. Perform a function test.

Figure 29: Adjusting the vertical lift on the ACROBAT LCH spring arm

9.4.3 Adjusting the vertical lift on the ACROBAT LCH spring arm

(See "Figure 29")

Tool to be used

Use the pin (Ø 4 mm x 110 mm) available as an accessory included in the scope of delivery.

Dismantling the cover panels

1. Unscrew 2 countersunk Phillips screws M3 x 6 mm ②.
2. Gently disengage each cover panel ⑤ on the 2 staggered straps ④ using a suitable slotted screwdriver.
3. Remove the 2 cover panels ⑤.

Adjusting the vertical lift

4. Insert the pin (Ø4 mm x 110 mm) into the perforated adjusting nut in the adjusting aperture ①.

To reduce the swivel range

- Turn the adjusting nut ① upwards in the anti-clockwise (-) direction as illustrated in the Figure.


To extend the swivel range

- Turn the adjusting nut ① downwards in the clockwise (+) direction as illustrated in the Figure.


Mounting the cover panels

5. Mount the cover panels ⑤ as described in "Chapter 9.3.3" on page 99.
6. Perform a function test.


10.1 ACROBAT DB 2001 ceiling-mounted version, single-arm variant
ACROBAT LCH/NRH


ACROBAT 2000 (standard version)


10.2 ACROBAT DB 2002 ceiling-mounted version, dual-arm variant
ACROBAT LCH/NRH


ACROBAT 2000 (standard version)

10.3 ACROBAT wall-mounted version, single-arm variant
ACROBAT LCH/NRH

ACROBAT 2000 (standard version)

Modes of operation	The pendant system is suitable for continuous operation.
Rating plates	The rating plates are attached to: <ul style="list-style-type: none"> • the top side of the extension arm • the top side of the spring arm • the approved adaptations
Dead weights of the intermediate ceiling kit	Intermediate ceiling rack (ceiling panel, 3 tubes and 3 threaded rods).....58.2kg Interface plate30.0kg
Dead weight of the flange with ceiling tube, variable ceiling tube length	Ceiling tube length 200 / 250 / 300 / 400 / 500mm..... 11.4 / 12.0 / 12.5 / 13.6 / 14.7kg Ceiling tube length 600 / 700 / 800 / 900mm 15.8 / 16.9 / 18.0 / 19.1kg Ceiling tube length 1000 / 1100mm.....20.2 / 21.3kg Canopy4.0kg
Dead weights of the ACROBAT 2000 ceiling-mounted version	Ceiling mount fixture – single-part (DB 2001) 10.6kg Ceiling mount fixture – dual-part (DB 2002)..... 13.9kg Extension arm 600 / 750 / 950mm 2.6 / 3.2 / 4.2kg Spring arm (standard version)6.2kg Spring arm NRH5.5kg Spring arm LCH5.7kg
Dead weights of the ACROBAT 2000 wall-mounted version	Wall bearing 1.8kg Extension arm 300 / 600 / 750/950mm 1.3 / 2.6 / 3.2 / 4.2kg Spring arm (standard version)6.2kg Spring arm NRH5.5kg Spring arm LCH5.7kg
Dead weight of the adaptations	Connecting piece 1kg Bracket.....5kg VESA adapter2kg See <i>“Chapter 12” on page 107.</i>
Maximum loading capacity on the ACROBAT 2000 extension arm	ACROBAT 2000 ceiling device 600mm extension arm21.0kg 750mm extension arm 18.0kg 950mm extension arm 15.0kg ACROBAT 2000 wall-mounted device 300mm extension arm21.0kg 600mm extension arm 15.0kg 750mm extension arm 13.0kg 950mm extension arm 10.0kg
Maximum loading capacity on the spring arms	In accordance with the spring arm equipment
Possible spring equipment in the spring arm (standard version)	1.5 - 3.5; 3.5 - 7; 7 - 12; 9 - 15; 12 - 18; 13.5 - 21 kg (see rating plate)
Possible spring equipment in the NRH spring arm	95 - 160; 95 - 176Nm (see rating plate)
Possible spring equipment in the LCH spring arm	95 - 160; 95 - 176Nm (see rating plate)
Electrical data	Depending on the customer-specific equipment (see rating plate)

Optional brakes	Friction brakes on the extension arm
Brakes	Friction brakes on the NRH and LCH spring arms
Noise level	Sound energy > 60db(A) (EN ISO 3744) are not exceeded
Operation	Manual forces < 550N
Medical Device Directive	Classification I
Protection class / type	ACROBAT 2000: Protection class in accordance with EN 60601-1 I ACROBAT LCH: Protection class in accordance with EN 60601-1 I IP classification in accordance with IEC 60529 IP 20
Applicable standards	<ul style="list-style-type: none"> • Medical Devices Act (Medizinproduktegesetz - MPG); • MDD 93/42/EEC, 2007 – Medical Device Directive; • EN 60601-1: 2006 - Medical Electrical Equipment - Part 1: General Requirements for Basic Safety and Essential Performance.
Approvals of the standard equipment	Recognised UL component.
Approved adaptations	<p>The following Ondal products are approved as adaptations to the pendant system:</p> <ul style="list-style-type: none"> • Ondal spring arms in accordance with “Chapter 11” on page 105 • Ondal adaptations in accordance with “Chapter 12” on page 107 • The components are adapted to each other and safe to operate. Any other type of installation, and in particular the use of components from third-party manufacturers, is strictly prohibited because these components can be potential sources of danger. • The combination of any other Ondal product with the pendant system must be approved by Ondal Medical Systems GmbH. If applicable, the conformity evaluation must be renewed.
Read the operating instructions for combined medical products.	<ul style="list-style-type: none"> • The pendant system can be equipped with adaptations and end devices of third-party manufacturers. To prevent dangerous overload, which can damage or lead to a collapse of the pendant system, the maximum loading capacities specified in “Chapter 11, Technical Data, on page 105” must be adhered to. • The party placing the appliance into operation is responsible for the validation of the overall system. A conformity evaluation procedure shall be executed if required and a declaration in accordance with Article 12 of 93/42/EEC (Medical Device Directive, MDD) shall be provided. • Read the Operating Instructions provided by the third-party manufacturer and in particular the relevant pages with information on the operation of the end device.
CE marking	Ondal Medical Systems GmbH • Wellastraße 6 • D-36088 Hünfeld • Germany

Designation	Dead weight	Maximum payload
CEMOR A1-N monitor connecting piece for mounting a flat screen	22kg	10kg
CEMOR A1-L monitor connecting piece for mounting a flat screen	3kg	10kg
CEMOR A1-F monitor connecting piece (pre-assembled) for mounting a flat screen	0.6kg	7kg
Universal connecting piece (pre-assembled)	0.3kg	21kg
Vario View flat screen support, single connection	4kg	12kg
Vario View flat screen support, dual connection	8.5kg	12,5kg
To order accessories please contact our sales team.		

Inspection Plan for

- Acrobat Swing Acrobat 2000 Acrobat 3000
- Acrobat 77 Acrobat LCH OndaSpace

System Data

Supplier: _____ Date of installation: _____
 _____ Ondal serial number: _____
 _____ Operator serial number: _____
 _____ Device location: _____

Important Information

- Inspections must be performed by trained service personnel.
- The inspection intervals must be observed.
- This inspection plan is only valid when combined with the Ondal Installation and Operating Instructions which must be complied with complementary to the inspections.
- After 10 years, the functional inspection must be performed once a year.

At the specified intervals, the pendant system must be inspected for the following points by Ondal, a company authorised by Ondal, or by personnel with the corresponding qualification:

Visual inspection (once a year)

- The parts of the pendant system are not deformed*
- The system is free from defects in paint wor
- The plastic parts are available and in position*
- The plastic parts are free from cracks*
- All rating plates are in position and legible

Period of use in years									
1	2	3	4	5	6	7	8	9	10
OK NOK	OK NOK	OK NOK	OK NOK	OK NOK	OK NOK	OK NOK	OK NOK	OK NOK	OK NOK
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Functional inspection (every 2 years)

- Free rotatability/limit stops (depending on the version)**
- horiz. a. vert. joints smooth running, grease if required
- Vertical lift correct, readjust if required**
- Check and grease securing segment*
- Securing ring in position and shape (extension/spring arr
- Load comp./spring tension correct, readjust if required
- Collision damage
 - all welding points free from cracks**
- Inspection of protective conductor transition resistance** (only applicable if current conducting cables are integrated)
- After 10 years, the cables and hoses (if existing) must be inspected and replaced if required**

2	4	6	8	10
OK NOK	OK NOK	OK NOK	OK NOK	OK NOK
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Confirmation of the inspection performed

The inspection tasks mentioned above were executed incl. the necessary adjustments and safety inspection:

1st year _____ Date 2nd year _____ Date 3rd year _____ Date 4th year _____ Date 5th year _____ Date	Signature/stamp _____ Signature/stamp _____ Signature/stamp _____ Signature/stamp _____ Signature/stamp	6th year _____ Date 7th year _____ Date 8th year _____ Date 9th year _____ Date 10th year _____ Date	Signature/stamp _____ Signature/stamp _____ Signature/stamp _____ Signature/stamp _____ Signature/stamp
---	---	--	---

* Damaged or deformed components should be replaced as a precaution. Plase contact the spring arm supplier.
 ** If one of these problems is detected during the inspection, the spring arm must be taken out of service immediately as a precaution in order to prevent any further damage to persons or equipment. Inform the system suppliers immediately.

The medical device handbook which belongs to each medical device and is prescribed in accordance with the German Medical Product User Regulations (MPBetreibV) must be available on site. Service and maintenance work and also safety inspections shall be documented in this medical device handbook. Inspection reports, including the present report, must be filed in the corresponding medical device handbook.

1540866, Edition 2012-02, Version 2

Subject to change and errors. Edition 2012-02